

MANDINGO® CONDOMS


Measurement Ruler

Please measure your penis properly by following instruction below:

To ensure that your Measurement Ruler printed correctly, and avoid any uncertainty when using the printed Measurement ruler, please refer to step 1. In case of printing issues, disable printer page scaling or refer to www.mandingocondoms.com FAQ's page for further information.

STEP 1

To ensure your ruler has printed to the correct size, placing a bank card onto the illustration below. If the card fits then page scaling is correct.



STEP 2

Please cut out the ruler.

STEP 3

To determine the length of your penis, measure the top side of an erect penis, not the side of your testicles. **Hold the ruler against the point where the shaft meets the pubic area (the base of your penis),** and simply bend the ruler along the shaft to the tip of the gland (**head**), and read measurement.

STEP 4

To determine the girth of your penis, you wrap measurement strip one inch up from the base of an erect penis to find your suggested size.



In case there is not a recommended condom size for you, please [contact us](#) for additional assistance.

Please note: The size determined for you is a suggested size. You may want to choose a size (One size larger or one size smaller) based on your individual preference.

ATTENTION: Please pay attention to sharp cutting edges during the application! User assumes all risk of bodily injury associated with usage.

